

Fona Opt digital FONA XPan 3D TAC Radiologische Röntgenstrahl



Valutazione: Nessuna valutazione

Prezzo:

Prezzo scontato: 37000,00 €

Prezzo di vendita scontato:

Prezzo di vendita: 37000,00 €


[Fai una domanda su questo prodotto](#)

Produttore: [FONA](#)

Descrizione



With XPan 3D you can see the smallest details for precise diagnostics. 3D voxel size of only 160 µm combined with 14 bit dynamic range, which equals more than 16.000 gray-levels, ensures high image quality after each exposure. See what you could not see before with XPan 3D.



With just one exposure you get the full arch field of view 8,5 x 8,5 cm. The reconstruction takes only 10 seconds, thanks to modern 64 bit technology, saving your time and the time of your patients. See the full dentition, mandibular canal and maxillary sinus in one single exposure and in record time.



TECHNICAL SPECIFICATIONS

Model	GM02 Flat Panel
Resolution	1024 x 1024
Field of View	8,5 x 8,5 cm
Reconstruction Time	10 seconds
Bit Depth	16 bit
Image Quality	High
Integration	Compatible with all major dental X-ray systems
Software	Proprietary software for dental X-ray reconstruction
Manufacturer	Autoklav Steri